



Eye Piece	Location where you view your specimen
Objective	Contains the different lens powers
Stage Clip	Used to hold slide on stage
Stage	Location to set the slide
Diaphragm	Controls the amount of light exposed to the slide
Rheostat	Controls the intensity of light
ON-OFF Switch	Turns light on or off
Course Adjustment	Moves stage closer or further from the objective lens quickly
Fine Adjustment	Moves stage closer or further from the objective lens slowly
Stage Clip Up-Down	Moves slide on stage up or down
Stage Clip Right Left	Moves slide on stage right or left
Slides to View	Letter e, onion root tip, bacterial capsule, cheek smear
Instructions	Location to learn how to use virtual microscope
Getting Started	Tutorial for learning how to use the virtual microscope
Checklist	Document to keep track of what you have learned
Switch Views	Allows user to view specimen in the microscope